

"Every breath you take"

N5 Biology, Unit 2, KA5: Transport systems - Plants.

Aim: To compare the stomatal density on the upper and lower surface of Tradescantia leaves.

Tradescantia provide a beautiful specimen from which to observe stomata: the deep purple colour of the leaf contrasts to the green stomata for clear observation. Numeracy and investigative skills are important in this practical activity.






Materials

- Light microscope
- Thick acetate
- Graph paper (mm squares)
- Photocopier
- Scalpel
- Tradescantia zebrina leaf.



Stage 1: Preparation of the graduated "slide"

Photocopy mmsquared graph paper onto thick acetate.






Using a microscope, observe the graduated "slide" and determine the area of the field of view at each magnification.




Cut the acetate into microscope slide sizes to make graduated "slides".



Stage 2: Observe stomata on the lower surface of a Tradescantia leaf


Using a scalpel, remove a small section of a leaf from Tradescantia.





Place the leaf on a microscope slide with the lower surface facing the objective lens

- Observe the stomata and count the number of stomata you can observe in the field of view. Note the magnification used to perform the count.
 - Calculate the number of stomata per squared centimetre of the Tradescantia.



x40 magnification



x100 magnification

How to determine the density of stomata in Step 4

- At x100 magnification, there were 13 stomata in the field of view.
- Using the "graduated slide", the diameter of the field of view at x40 magnification was 4 mm.
- Diameter of the field of view at x100 magnification:
 - \circ (40/100) x 4 = 1.6 mm
- Area of the field of view at x100 magnification:
 - $\sigma = \pi r^2 = \pi x (0.8)^2 = 2.01 \text{ mm}^2$
 - o 13 stomata / 2 mm²
 - o 6.5 stomata / mm² of lower leaf surface

Stage 3: Observe stomata on the upper surface of a Tradescantia leaf

Repeat Stage 2 for the upper surface of the leaf.



Calculate the percentage change between the upper and lower surface of the leaf.



Questions to explore

Suggest a plant that might not show this distribution in stomata.

Explain the conclusion you have drawn from your data.