

PROTOCOL INVESTIGATING PHOSPHATASE ACTIVITY



sserc

*Extract phosphatase
from bean sprouts and
catalyse a
degradation reaction*

Page

1

Background

2

Risk Assessment

3

Materials

3

Method

5

Results

6

Supplementary resources



Background

Phosphatase enzymes are involved in a range of metabolic reactions. A key function of these enzymes is to release phosphate groups into the metabolic pool thereby increasing their availability for use in a range of processes including ATP synthesis and membrane construction.

Acid phosphatases (those with an optimum pH <7.0) can be extracted from a range of plant tissues - germinating mung beans or bean sprouts are a cheap and reliable source. The substrate is phenolphthalein bisphosphate (PPP). Under suitable conditions, phosphatase catalyses the breakdown of PPP to form phenolphthalein (PP) (Figure 1).

At neutral or acidic pH, the products of this reaction (PP and phosphate) are both colourless - so their presence is difficult to detect. This can be overcome through the addition of sodium carbonate which has 2 effects:


1) raising the pH of the solution >pH10 with cessation of enzyme activity,

2) converts PP to its anionic form, which is pink.



*NMR structure of a
phosphatase enzyme from
Arabidopsis thaliana*

Figure 1: Phosphatase-catalysed breakdown of PPP to form free phosphate and PP, both of which are colourless. The addition of sodium carbonate stops further enzyme activity and converts PP to its anionic, pink form.



AIM

To investigate the effect of pH on phosphatase activity in bean sprouts.

This protocol can be adapted to investigate a range of independent variables, including:

- effect of temperature on phosphatase activity
- effect of enzyme concentration on phosphatase activity
- effect of substrate concentration on phosphatase activity
- effect of end-product inhibition on phosphatase activity
- effect of tissue type on phosphatase activity.

RISK ASSESSMENT

A risk assessment for this activity can be downloaded from the SSERC website. Click [here](#). This should be adapted for your centre, where appropriate.

Briefly, the main hazard associated with this protocol is the use of a centrifuge. This should be PAT tested and care must be taken to ensure the lid cannot be opened while the rotor is spinning. The centrifuge tubes must be accurately balanced in the rotor.



MATERIALS REQUIRED PER PAIR

Part 1 - Preparation of the enzyme extract

- 20 g bean sprouts
- pestle and mortar
- water
- 3 cm³ plastic pipettes
- microfuge
- 6x microfuge tubes
- marker pen
- container to store enzyme extract

Part 2 - Phosphatase Assay

- enzyme extract
- 25 cm³ 10% (w/v) sodium carbonate
- stopwatch
- access to a waterbath (30°C)
- paper towels
- 14x absorption cuvettes
- colorimeter (550 nm)
- cuvette rack
- 10 cm³ citric acid / phosphate buffer (pH 5.0)
- 10 cm³ citric acid / phosphate buffer (pH 7.0)
- 6 cm³ 0.2% phenolphthalein phosphate
- 1 cm³ automatic pipette and tips
- polystyrene cup

OVERVIEW OF METHOD

1. Crush 20 g bean sprouts in a mortar with 5 cm³ water.



2. Divide extract between 6 microfuge tubes.



3. Centrifuge for 5 minutes. Store supernatant for next step.



4. Incubate buffer, enzyme and substrate at 30 °C.



5. Add 1 cm³ sodium carbonate to 7 cuvettes.



6. Mix 2 cm³ enzyme, 10 cm³ buffer, 2 cm³ substrate. Every 2 minutes, transfer 1 cm³ to a cuvette



STEP-BY-STEP METHOD

Part 1 - Preparation of the enzyme extract


1. Place 20 g bean sprouts in a mortar. Remove and discard the green testa (seed case) if it is attached. Germinated mung beans can also be used - use about 30 mung beans per extraction.
2. Add 5 cm³ water to the mortar. Grind the bean sprouts with a pestle to achieve a smooth paste.
3. Cut the tip off a plastic pipette. Divide the extract equally between 6 microfuge tubes (should be approximately equal volume to ensure the centrifuge rotor is balanced).
4. Centrifuge the samples for 5 minutes. Using a plastic pipette, transfer the supernatant from the microfuge tubes to a labelled container.

To carry out the assay at pH 5 and pH 7, at least 5 cm³ enzyme extract is required.



Part 2 - Phosphatase Assay

The following diagram gives an overview of this part of the protocol.



All solutions should be kept at 30 °C throughout the assay. It is accurate to dispense the volumes outlined in the following steps using an automatic 1 cm³ pipette with a clean tip.

- 1 Transfer water from the waterbath (30 °C) to the polystyrene cup. Stand the container of buffer (pH 5), enzyme extract and substrate (PPP) in the polystyrene cup. This provides each learner with their own individual waterbath. If carrying this out for an assessed piece of work, a thermostatically-controlled waterbath should be used instead.

In the image opposite, two polystyrene cups have been used, with water at 30 °C. The substrate and enzyme are stored in bijou bottles in one cup (left), while the buffer is held in a universal bottle in the second cup (right).




- 2 Using the 1 cm³ automatic pipette, add 1 cm³ sodium carbonate into 7 cuvettes. Discard the tip.

- 3 Transfer 2 cm³ enzyme extract to 10 cm³ buffer solution.




- 4 For the colorimetric blank: transfer 1 cm³ of the mixture from the container of buffer + enzyme (step 3) to the first cuvette (containing sodium carbonate).



5 Transfer 2 cm³ PPP to the enzyme/buffer mixture. Start the stopwatch. Mix the contents thoroughly but without creating too many bubbles.



6 At 2 minute intervals, transfer 1 cm³ of the PPP/ enzyme / buffer mixture to a cuvette. The sodium carbonate will stop the reaction and convert the product (PP) to its anionic form, which is pink.



7 Use the colorimetric blank (step 4) to zero the colorimeter. This contains buffer and enzyme only and is held in cuvette 1.

Measure the absorbance of the remaining solutions at 550 nm. If using the Mystrica colorimeter, this will require using the green diode.

The assay can now be repeated using the buffer at pH 7.0.



RESULTS

Present results in a table similar to the one opposite.

The results can be then plotted as a line graph.


Sample results are shown on the following page.

Time (min)	Absorbance	
	pH 5.0	pH 7.0
0		
2		
4		
6		
8		
10		
12		




Time (min)	Absorbance	
	pH 5.0	pH 7.0
0	0.00	0.00
2	0.09	0.02
4	0.18	0.04
6	0.22	0.07
8	0.31	0.12
10	0.42	0.19
12	0.56	0.26


Effect of pH on phosphatase activity in bean sprouts



When this experiment is carried out over a range of pH values, the following activity is observed (right). This experiment involved taking a fixed time-point measurement after 10 minutes at each pH value, showing an alternative approach to this protocol.



This chart (right) presents expected results from the protocol when investigating temperature as the independent variable.



SUPPLEMENTARY RESOURCES

Fun with Phosphatase

In a previous issue of the SSERC Bulletin (1) we explored how the enzyme *alga oxidase* might be used to support learning and teaching across a range of curriculum levels. Phosphatase is one of the enzyme systems specifically mentioned in both of the recently implemented CfE Highers in Biology (2) and Human Biology (3). A particular emphasis is the use of phosphatase in experiments on product inhibition.

Phosphatase enzymes are involved in a range of metabolic reactions. A key function of these enzymes is to release phosphate groups into the metabolic pool, thereby increasing their availability for use in a range of processes including ATP synthesis and membrane transport. Acid phosphatases (those with a preference for acidic pH) can be active over a range of pH values with particularly strong activity being particularly characteristic of those from the acidic environment of lysosomes.

Under suitable conditions the enzyme phosphatase will catalyse the following reaction in these substrates:

An acidic pH will shift the equilibrium of the above reaction towards the products (P_i and phosphate) and this process is reversible, allowing the reaction to be reversed.

Insoluble through the addition of sodium carbonate which has a pH effect to raise the pH of the solution to above 10 with consequent cessation of enzyme activity and CO₂ evolution. Insoluble precipitates will be formed from both the substrate and the product in solution.

SSERC

Phosphatase Assay – Technical/Tutor Guide

1.0 Background

The experiments described here are based upon information found on the Science and Plants for Schools (SAPS) website (see www.saps.org.uk) and a publication by Barry Meatyard (Phosphatase enzymes from plants. *Journal of Biological Education*, 33 (2), 109-112).

Technician Guide

SESSION 1

INVESTIGATING PHOSPHATASE ACTIVITY

[PPT to download for lessons](#)

SSERC Risk Assessment (revised version March 2018)
(based on HSE's INDG 163 'Risk assessment - A brief guide to controlling risks in the workplace')

2 Pitreavie Court, South Pitreavie Business Park, Dunfermline KY11 8JU
tel: 01383 636076 e-mail: gs@sserc.co.uk web: www.sserc.org.uk

Activity/assessment	Extraction of phosphatase from bean sprouts and subsequent degradation reaction		
Date of assessment	25th July 2023		
Date of review (Step 5)	25th July 2024		
School	SSERC		
Department	Biology		

Step 1	Step 2	Step 3	Step 4
List Significant hazards here:	Who might be harmed and how?	What are you already doing? What further action is needed?	Actions by whom? Due date Done
Centrifuge	Technician, Teacher or Learner	Centrifuge must be suitably maintained with a valid PAT test certificate.	
Hazard – injury to operator by spinning rotor or electrical fault.		The design should include an interlock to prevent the lid being opened while the rotor is turning.	
Hazard – electrical fault		The centrifuge tubes must be accurately balanced and positioned opposite each other in the rotor.	

SSERC bulletin (2015) available to [download](#).

[Risk Assessment](#)